

Pro-Wrap Application

PW110 - PW190

- Fuss-free from ice cubes, bandages and hot water bottles
- Reusable Cold and Hot applications
- Breathable and ergonomic brace
- Uniquely designed for each specific body
- Unique combination for uses in therapy clinics, rehabilitation centers, athletic facilities at home
- 30 minutes auto shut off safety device

PW110 Shoulder and Upper Back

- Great for both acute and chronic shoulder injuries
- Fits either shoulder
- Suitable for chest size within 50"



PW120 Elbow

- Great for both acute and chronic elbow injuries



PW130 Wrist

- Great for both acute and chronic wrist injuries
- Fits either hand



PW140 Lower Back

- Great for both acute and chronic lower back injuries
- Suitable for waist size within 46"



PW150 Thigh

- Great for both acute and chronic thigh injuries



PW160 Knee

- Great for both acute and chronic knee injuries



PW170 Ankle

- Great for both acute and chronic ankle injuries



PW180 Foot

- Great for both acute and chronic ankle injuries
- Fits either foot



PW190 Shoulder and Upper Arm

- Great for both acute and chronic shoulder injuries
- Fits either shoulder
- Suitable for chest size within 50"




Product Specification

Series	ProWrap				
Model	PW110	PW120	PW130	PW140	PW 150
Parts	Shoulder	Elbow	Wrist	Lower Back	Thigh
Size	< 50"	one size	one size	< 46"	one size
Dimension (mm)	860x590	570x260	500x95	1100x200	690x200
Power (W)	14	14	7.5	16	16
Cabric® Heating Pad Size (mm)	250x150	250x150	260x70	320x160	320x160
Neo Cold Pack size (mm)	250x150	250x150	260x70	310x160	310x160
Series	Pro Wrap				
Model	PW160	PW170	PW180	PW190	
Parts	Knee	Ankle	Foot	Shoulder	
Size	one size	one size	one size	< 50"	
Dimension (mm)	630x270	260x290	320x220	860x590	
Power (W)	11	12	7.5	14	
Cabric® Heating Pad Size (mm)	240x220	220x210	260x70	280x240	
Neo Cold Pack size (mm)	230x220	215x215	260x70	280x240	
Temperature Controller	> One-Button operation > 3 settings temperature control > 30 minutes auto shut-off				
Adaptor	> Input: 100-240V / 50-60HZ > Output: 12VD / 1.5A > medical standard certified				
Brace / Wrap	> Neoprene > Nylon > Other				
Neo Cold Pack	> Exterior: PVC > Fillin Liquid: Nontoxic gel				
Cabric® Heating Pad	> Exterior (Black): SBR / Nylon > Heating Fabric: 100% Carbon Fiber Fabric > Lining: Tricot Fabric				
Car Charger	> Input / Output: 12V only (optional accessory)				
Color	Orange / Black				

FL_PW_1_4011

Pro-Wrap Series

3 in 1 HOT · COLD · BRACE Therapy



HOT Therapeutic Far Infrared Energy reduces muscle tension and joint stiffness

COLD FIR energy travels in straight lines and is easily absorbed by the human body

BRACE Ergonomic and breathable braces provide more support and comfort on your 8 main body parts

CONTROLLER 3 setting temperature with auto power off time precise temperature control.



TherMedic Pain Management

Healthy living is based on everyday care and prevention. Our research team has found that current pain management programs tend to focus on treatment while ignoring prevention. The importance of prevention is often only appreciated after pain strikes. Our main goal at TherMedic is to increase self-awareness regarding pain and to create a pain-free quality of living. By building on a foundation of preventative medical care, we integrate people-oriented mentalities into our products to provide new methods of household healthcare. All of our products are exclusively equipped with cabric heating pads, made of aerospace grade carbon fiber fabric, which allow far-infrared waves to effectively and efficiently reach deep tissues via resonance absorption. This boosts circulation, which is the cornerstone of self-healing. Through Thermedic, heat and health come as easily and naturally as the air we breathe, and we hope the whole world will be able to feel our determination and enthusiasm.

Pain is most classified as being either acute or chronic. Pain management is the medical discipline concerned with the relief of pain. Conquering all the conventional inconveniences, TherMedic is born for the purpose of improving the quality of human life. With benefits of pain relief, improvement of blood circulation and promoting healing, TherMedic saves you from pain killers and expensive physiotherapy visits and takes care of your bothering pain anywhere at home, office and outdoor with total solution of pain management.

Far Infrared Technology



More than 15 years of experience in the field of carbon materials have led TherMedic to the development of cabric, a carbon fiber fabric with the ability to emit FIR (far infrared) energy. Cabric's state of the art production facilities in combination with highest quality standards make this fabric a unique base material for a variety of different applications, including medical therapy.


Compared to conventional heating pads which only produce "line heat", the TherMedic far infrared product line is using the superior far infrared technology (FIR). Far infrared is a scientifically proven form of heat therapy that works more efficiently and effectively than traditional heating techniques. Most notably, far infrared is a completely safe form of thermal radiation, as it does not emit the harmful electromagnetic frequencies (EMF) like most other electronic products.

Unlike conventional heating technology which is using metal wire as heating elements, the unique cabric inside technology, a carbon fiber fabric which in itself is electrically conductive and therefore does not use any wires. The fabric naturally emits far infrared energy which supports various therapeutic applications. No matter if you are suffering from acute or chronic pain or just use it for preventive therapy, TherMedic's extensive product line is the right choice.

3 in 1 THERAPY


3in1 THERAPY - HOT · COLD · BRACE

In this ever-changing and high competitive world, everyone desires a pain-free. Pro-Wrap plays a key role as your personal pain manager providing a 3 in 1 therapeutic solution.




HOT THERAPY

Made of 100% aerospace grade carbon fiber fabric, flexible cabric heating pad with benefits of even heating pad therapeutic Far Infrared Energy reduces muscle tension and joint stiffness, penetrates to deeper tissues to boost blood circulation, increase metabolism and accelerate healing.



COLD THERAPY

Neo Cold Pack narrows vessels and decreases blood flow, reducing swelling and inflammation within the first 24 - 48 hours of the acute stage of an injury to prevent tissue damage.




BRACE

Ergonomic and breathable braces provide more support and comfort on your 8 main body parts with 100% coverage protecting from second injury.

Features

- Great for both acute and chronic injuries.
- Reusable Cold and Hot applications.
- Breathable and ergonomic brace.
- Uniquely designed for each specific body part.
- Unique combination for uses in therapy clinics, rehabilitation centers, athletic facilities and at home.



How to use?

Set up for HOT THERAPY


- Insert Heating Pad into the mesh pouch
- Connect Heating Pad connector with temperature controller
- Plug in AC / DC Adaptor
- Push the „Power on“ button

Set up for COLD THERAPY

- Take out Heating Pad from the mesh pouch
- Pull out the connector carefully
- Insert Neo Cold Pack into the mesh pouch

Set up for BRACE

- Remove Heating Pad and Cold Pack from the mesh pouch, simply use the Wrap as a brace.



1-9 Shoulder

- Slap Tear
- Dislocated Shoulder
- Labral Tear
- Rotator Cuff Bursitis
- General Shoulder Pain
- Frozen Shoulder

2 Elbow

- Golfer's Elbow
- Medical Epicondylites
- Tennis Elbow
- Lateral Epicondylites
- Bicep Tendon Rupture
- General Elbow Pain

3 Wrist

- Carpal Tunnel Syndrome
- DeQuervain's Tenosynovitis
- Wrist Sprains

4 Lower Back

- Sprained Back Muscles
- Sciatica
- Muscle Spasms
- Lumber Pain

5 Thigh

- Hamstring Pulls
- Muscle Strains
- Surgery Recovery
- Upper Quadriceps Strains

6 Knee

- Mensical Injuries
- Osteoarthritis
- Bursitis in the Knee
- Osgood-Schlatter Disease (OSD)
- Patellar Tendonitis
- Anterior Cruciate Ligament (ACL) Tear
- Posterior Cruciate Ligament (PCL) Tear
- Lateral (LCL) & Medial Collateral Ligament (MCL) Tear

7 Ankle

- Achilles Tendon Tear
- Achilles Tendon Inflammation
- Achilles Tendon Tenderness
- Ankle Sprains
- Ankle Strains
- Twisted Ankle

8 Foot

- Plantar Fasciitis
- Heel Spurs
- Tendonitis
- Tendinopathy in the foot
- Arch Pain
- Turf Toe