


General Catalogue

Kenz



SUZUKEN COMPANY

SINGLE CHANNEL DIGITAL ELECTROCARDIOGRAPH

KENZ ECG 108/110



Kenz ECG 108

COMPACT AND LIGHTWEIGHT

Less than 2.5 kg including the built-in rechargeable battery

EXCELLENT RECORDING QUALITY

Clear and precise ECG tracings

Digital filters eliminate AC interference, baseline drift and muscle tremor (EMG).

EASY OPERATION

Recording by just one touch of the RUN/STOP key

ALPHANUMERIC PATIENT INFORMATION ENTRY

Patient ID number (up to 12 digits), sex, name and comments (up to 36 letters) can be entered by using the operation keypads.

And an external PS/2 keyboard can be connected to enter this information.

RECHARGEABLE BATTERY OR AC OPERATION

Maximum one hour continuous operation with fully charged for 2.5 hours (Battery operation)

ECG PREVIEW DISPLAY (ECG 110 ONLY)

Large back-lit LCD (5.7 inch)

MESSAGE DISPLAY

VARSATILE RECORDING MODES



Kenz ECG 110

3 CHANNEL DIGITAL ELECTROCARDIOGRAPH WITH ANALYSIS

KENZ CARDICO 302

ECG preview display on 5.7 inch back-lit LCD

Display freeze and print copy

Highly accurate computerized interpretation based on simultaneous acquisition of all 12 leads

Rechargeable battery and AC operation

Data transmission to ECG data management software, ECG Data Viewer CVS-02

Versatile 6 modes of operation

- Auto mode
- Post-stress test mode
- Freeze mode
- Manual mode
- Pre-check mode
- Long-term mode



Kenz Cardico 302

12 CHANNEL DIGITAL ELECTROCARDIOGRAPH WITH ANALYSIS

KENZ CARDICO 1210



- | 3, 4, 6 or 12-ch ECG recording with analysis
- | Auto, Manual and Long-term recording
- | 1 sheet, 2 sheets or 3 sheets printout
- | Compartment for A4 size Z-fold paper
- | 5.7 inch LCD display
- | Soft keypad (option:external keyboard)
- | Built-in memory for 100 ECG data
- | Power by AC or DC
- | Excellent recording quality
- | Data transmission to Kenz ECG data viewer CVS-02

Network ECG System Cardico 1211

ECG management encompassing data management and network environment are now required in Physician's Offices, Clinics and Hospitals more than ever.

| High Resolution, Wide Viewing Angle Touch Screen Colour Display

| Large Memory

1000 records in internal memory and additional 1000 records can be stored in the optional filing card.

| Network ECG

Cardico 1211 comes with a Local Area Network (LAN) port that enables functional patient data management and electronic filing with the optional Kenz ECG Data Viewer CVS-02.

Wireless Transmission is possible by using locally purchased LAN module.

| New Arrhythmia Screening

The unique Long-term function features the full disclosure and summary report for Supra-Ventricular Premature Contraction. Ventricular Premature Contraction and Atrial Fibrillation.

| The multi language version

Chinese, English, French, Germany, Italian, Russian and Spanish




ECG Data Management System

Kenz ECG Data Viewer CVS-02

Powerful & Versatile Software that Makes ECG Data Management simple as possible

All of Kenz multi-channel ECG, Cardico series, has connectivity for paperless data management system through the ECG Data Viewer CVS-02.

- Edit interpretation statement and comments
- Arrhythmia analysis up to 5 minutes and ANST (Autonomic Nervous System Test) report
- Import other vital sign data from other devices and save as PDF
- 12 lead ECG display, over-read and serial ECG comparison of patient's data
- Save in graphic image format [JPG, PNG, BMP, PDF]
- Import and export patient information




Analysis Report (Over-read 3x4)




12 lead ECG display



Arrhythmia Analysis display



Arrhythmia Analysis Report



ECG Network System

Digital Holter ECG & BP Analysis System

Super Speedy Scanning and Full Automatic Analysis

Kenz Holter ECG and ambulatory blood pressure system makes Holter diagnosis quick and easy for even General Practice. Compact size and water proof features of Cardy recorder improve patient's amenity ever.



Cardy 301



Cardy 302*



Cardy 301L*

Digital Holter ECG Recorder **Kenz Cardy 301**

- Compact recorder small as business card size and light weight 72 grams
- 24-hour recording with 2ch/3ch
- Durable solid state casing without LCD
- Full time pacemaker detection

*Cardy 301L with 48 hours recording and Cardy 302 water proof with disposable electrode will be launching soon.

Holter ECG Analysis Software **Kenz Cardy Analyzer 03**

- Speedy scanning and fully automatic analysis in 60 seconds*
*24-hour data
- Full automatic process allows simultaneous printing of Full Disclosure while downloading 24-hour data from Cardy 301 and archiving in about 5 minutes.
- Simple mouse-action editing for the final report AF detection and HRV display
- Connectivity with ambulatory blood pressure monitor BPM AM300CE



Kenz Ambulatory BPM AM300CE

Holter ECG Viewer and Card Writer **Kenz Cardy Controller 02**

- Patient demographic information transmission to Cardy 301/301L/302 using IrDA
- ECG preview on screen and remote start
- ECG searching by event, time for review



Kenz Cardy Controller 02

Blood and Infusion Warmer

Accurate & Reliable Temperature Management

Kenz Animec warmers are the solution to medical institutions for the most cost containment of warming disposables during operation and recovery stage of patient.



Kenz Animec AM301*

- | Electric plate dry heating for tube
- | No disposable bag required
- | Two Flow rate: 1-15 ml/min or 1-20 ml/min
- | Warming temperature setting: 36°C-39°C
- | Audible and visual temperature Alarm
- | Light weight and portable
- | Splash resistance

*Animec AM301 will be launching early 2009.

Kenz Animec AM4

- | Electric plate dry heating for disposable bag
- | Flow rate: 0-60ml / min (gravity flow)
0-160ml / min (with infusion pump)
- | Warming temperature: 31°C -37°C
- | Audible and visual temperature alarm
- | Splash resistance



Kenz Animec AM2S



- | Electric plate dry heating
- | No disposable bag required
- | Warming temperature:
27°C -37 °C
- | Flow rate: 1~12ml / min
- | Light weight and portable
- | Splash resistance

Kenz Flexwarmer mobile warmer



- | No electric source required
- | Reusable heating gel pack
- | No disposable bag required
- | Warming temperature:
27 °C -37 °C
- | Flow rate: 1~12ml / min
- | Super light weight
- | Operative heating time:
60 min/heating pack

Physical Activity Monitor

Not all steps are equal!

Lifecorder® technology uses a solid state medical-grade sensor and digital filtering to deliver extremely accurate results. Unlike an inexpensive pendulum sensor pedometer, the Lifecorder filters out motion and vibration and only counts real steps. It also measures the INTENSITY of physical activity and accurately records how much time is spent at moderate or greater intensity level activities and energy expenditure.

Lifecorder® e-STEP



Lifecorder e-STEP, the economical activity monitor for a general user, has the same internal mechanism and accuracy as other Lifecorder models, but without a computer download function. Like other Lifecorders, it measures step count, total calorie expenditure and keeps track of total activity minutes at moderate or greater levels of intensity in each 24 hour period with a 7 day total memory.



Lifecorder® PLUS



The new Lifecorder PLUS has a large graphic display which illustrates how step count, calorie expenditure, distance and activity minutes compare with personal daily targets. The bar graph shows the activity intensity level in 2 minute intervals over one hour periods. The Lifecorder PLUS is the ideal model for a user wanting to track real-time progress toward daily targets.



Lifecorder® EX



The compact Lifecorder EX, with a 200 day memory capability and battery life of six months, is especially designed for lifestyle management programs. It features two unique functions, including a key on the front face for logging exercise time, meal time, and medication time. The second exclusive function is the display of the real-time intensity value of physical activity for checking exercise levels in 4 second intervals.



Step Counts

Activity Minutes

Calories Expenditure

Activity Intensity Graph (METs)

Distance

Memory

Set Up

Target Set Up

Data Management

Kenz®

INNOVATIVE MEDICAL DEVICE & SOFTWARE



Kenz. is registered trade mark of Suzuken Company.

SUZUKEN COMPANY LIMITED

8 Higashikataha-machi
Higashi-ku, Nagoya
Aichi 461-0015 JAPAN
TEL: +81-52-950-6327 FAX: +81-52-962-7440
<http://www.suzuken-kenz.com>